



## IsomerFree™ RNA AMIDITES AND SUPPORTS: CONVENTIONAL (3' → 5' OLIGONUCLEOTIDE SYNTHESIS)



5'-DMT-2'-TBDMS-3'-Amidites

<b>B;</b>	Adenosine (n-bz)	ANP - 5671
	Cytidine (n-bz)	ANP - 5672
	Guanosine (n-ibu)	ANP - 5673
	Uridine	ANP - 5674
	Cytidine (n-ac)	ANP - 6676




5'-DMT-2'-TBDMS-3'-CPG

<b>B;</b>	Adenosine (n-bz)	N - 6101
	Cytidine (n-bz)	N - 6102
	Guanosine (n-ibu)	N - 6103
	Uridine	N - 6104
	Cytidine (n-ac)	N - 6106

## RNA SYNTHESIS IN REVERSE DIRECTION (5' → 3' -END)


- Coupling efficiency per step is very high (~99%).  
Results in Ultrapure RNA.
- Highly suitable for:  
*Attachment of 3'-Cholesterol using cholesterol amidites*  
*Attachment of 3'-PEGS using PEG amidites*
- Ideal for introducing lipophilic modifications & delivery reagents to the 3'-end of RNA oligonucleotides.
- A Number of applications are possible for easy attachment at the 3'-end of an oligonucleotide.
- Produced with an HPLC purity of greater than 98% and 31P NMR purity greater than 99%.
- Leads to smooth 3'-conjugation of macromolecules to Synthetic RNA.
- Reverse RNA synthesis results in complete absence of M + 1 species.

**Devoted to providing you with invaluable customer service and comprehensive technical support.**



3'-DMT-2'-TBDMS-5'-CPG

**B;** A(n-bz)  
C(n-bz)  
C(n-ac)  
G(n-ibu)  
U



3'-DMT-2'-TBDMS-5'-Amidites

B	Protection	Catalog #
A	n-bz	<b>ANP-3401 (Amidite)</b> <b>N-6201 (CPG)</b>
C	n-bz	<b>ANP-3402 (amidite)</b> <b>N-6202 (CPG)</b>
G	n-ibu	<b>ANP-3403 (amidite)</b> <b>N-6203 (CPG)</b>
U	n/a	<b>ANP-3404 (amidite)</b> <b>N-6204 (CPG)</b>
C	n-ac	<b>ANP-3405 (amidite)</b> <b>N-6204 (CPG)</b>

Request our new 2010 e-catalog  
[www.chemgenes.com](http://www.chemgenes.com)

Contact us for a quote today!  
**+1-978-694-4500**